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## STABILITY TEST REPORT

LYHER® Novel Coronavirus (2019-nCoV) IgM/IgG Antibody Combo Test Kit (Colloidal Gold)

### 1. Objective:

To study the stability of LYHER® Novel Coronavirus (2019-nCoV) IgM/IgG Antibody Combo Test Kit produced by Hangzhou Laihe Biotech Co.,Ltd.

### 2. Equipments:

- ☐ Incubators
- ☐ Thermometers
- ☐ Refrigerators

### 3. Materials

- ☐ LYHER® Novel Coronavirus (2019-nCoV) IgM/IgG Antibody Combo Test Kit finished products (hereinafter “LYHER® Kit”)
- ☐ Manufacturer’s reference products

### 4. Testing procedures

As described in the package insert of LYHER® Kit

### 5. Study plan

#### 5.1. Stability Test under Storage Condition (Real-time Stability)

Three batches of LYHER® Kit (140 pieces of each batch) are taken to store in QC laboratory. The temperature of storage condition is recorded every day. Products are ensured to be kept between 2°C and 30°C. If the temperature are to be lower than 2°C or higher than 30°C, the air conditioning/heating would be turned on to adjust the room temperature. Testings are scheduled to perform at 3, 6, 12, 18, 24, 30, 33 and 36 months. Furthermore, to test the stability of specimen diluents that are opened simultaneously with each batch, all specimen diluents remained open while sitting in room temperature for 30-45 minutes; then,

these specimen diluents are resealed and placed into an incubator.

Specimens of serum used for the stability test are listed below:

SN. Serum Description

N01~N05 Negative Both IgM and IgG are negative

P01~P05 Positive Both IgM and IgG are positive

M-L1~M-L5

Near IgM cut-off value

Prepared by diluting positive serum by negative serum according to a certain percentage. With the increase of dilution factor, the concentration

2

of positive specimens decreased gradually until can not be tested. We

took the last one IgM positive specimens that can be tested out (the color of the T line is G3) as the IgM cut-off value reference.

G-L1~G-L5

Near IgG cut-off value

Prepared by diluting positive serum by negative serum according to a certain percentage. While the increase of dilution factor, the concentration of positive specimens decreased gradually until can not be tested. We took the last one IgM positive specimens that can be tested out (the color of the T line is G3) as the IgG cut-off value reference.

Table 5.1-1: Results of Stability Test under storage condition (Month 0, tested on February 14, 2020)

Specimen

IgG Test Results IgM Test Results

	2002003	2002004	2002005	2002003	2002004	2002005
N01	Negative	Negative	Negative	Negative	Negative	Negative
N02	Negative	Negative	Negative	Negative	Negative	Negative
N03	Negative	Negative	Negative	Negative	Negative	Negative
N04	Negative	Negative	Negative	Negative	Negative	Negative
N05	Negative	Negative	Negative	Negative	Negative	Negative
P01	Positive	Positive	Positive	Positive	Positive	Positive
P02	Positive	Positive	Positive	Positive	Positive	Positive
P03	Positive	Positive	Positive	Positive	Positive	Positive
P04	Positive	Positive	Positive	Positive	Positive	Positive
P05	Positive	Positive	Positive	Positive	Positive	Positive
M-L1	NA	NA	NA	Positive	Positive	Positive
M-L2	NA	NA	NA	Positive	Positive	Positive
M-L3	NA	NA	NA	Positive	Positive	Positive
M-L4	NA	NA	NA	Positive	Positive	Positive
M-L5	NA	NA	NA	Positive	Positive	Positive
G-L1	Positive	Positive	Positive	NA	NA	NA
G-L2	Positive	Positive	Positive	NA	NA	NA
G-L3	Positive	Positive	Positive	NA	NA	NA
G-L4	Positive	Positive	Positive	NA	NA	NA
G-L5	Positive	Positive	Positive	NA	NA	NA

Table 5.1-2: Results of Stability Test under storage condition (Month 3, tested on May 14, 2020)

Specimen IgG Test Results IgM Test Results

	2002003	2002004	2002005	2002003	2002004	2002005
N01	Negative	Negative	Negative	Negative	Negative	Negative
N02	Negative	Negative	Negative	Negative	Negative	Negative
N03	Negative	Negative	Negative	Negative	Negative	Negative
N04	Negative	Negative	Negative	Negative	Negative	Negative
N05	Negative	Negative	Negative	Negative	Negative	Negative
P01	Positive	Positive	Positive	Positive	Positive	Positive
P02	Positive	Positive	Positive	Positive	Positive	Positive
P03	Positive	Positive	Positive	Positive	Positive	Positive
P04	Positive	Positive	Positive	Positive	Positive	Positive
P05	Positive	Positive	Positive	Positive	Positive	Positive
M-L1	NA	NA	NA	Positive	Positive	Positive
M-L2	NA	NA	NA	Positive	Positive	Positive
M-L3	NA	NA	NA	Positive	Positive	Positive
M-L4	NA	NA	NA	Positive	Positive	Positive

M-L5 NA NA NA Positive Positive Positive  
 G-L1 Positive Positive Positive NA NA NA  
 G-L2 Positive Positive Positive NA NA NA  
 G-L3 Positive Positive Positive NA NA NA  
 G-L4 Positive Positive Positive NA NA NA  
 G-L5 Positive Positive Positive NA NA NA

The real-time stability of Lyher kit was only tested for 0 month and 3 month, and the real-time stability will be continuously tracked.

### 5.2. Transport stability testing

Three batches of LYHER® Kit were mailed to Haikou City and Chifeng City; then, these kits were sent back to the manufacture facility. LYHER® Kits are randomly selected to test these specimens. Once the LYHER® Kit are expired, they will be randomly selected again to test specimens.

Table 5.2-1 Test results of the LYHER® Kits returned from Haikou City (Tested at March 6, 2020)

Specimen

Test results of LYHER® Kit when testing IgG Test results of LYHER® Kit when testing IgM  
 2002003 2002004 2002005 2002003 2002004 2002005

N01 Negative Negative Negative Negative Negative Negative

N02 Negative Negative Negative Negative Negative Negative

4

N03 Negative Negative Negative Negative Negative Negative

N04 Negative Negative Negative Negative Negative Negative

N05 Negative Negative Negative Negative Negative Negative

P01 Positive Positive Positive Positive Positive Positive

P02 Positive Positive Positive Positive Positive Positive

P03 Positive Positive Positive Positive Positive Positive

P04 Positive Positive Positive Positive Positive Positive

P05 Positive Positive Positive Positive Positive Positive

M-L1 NA NA NA Positive Positive Positive

M-L2 NA NA NA Positive Positive Positive

M-L3 NA NA NA Positive Positive Positive

M-L4 NA NA NA Positive Positive Positive

M-L5 NA NA NA Positive Positive Positive

G-L1 Positive Positive Positive NA NA NA

G-L2 Positive Positive Positive NA NA NA

G-L3 Positive Positive Positive NA NA NA

G-L4 Positive Positive Positive NA NA NA

G-L5 Positive Positive Positive NA NA NA

Table 5.2-2 Test results of the LYHER® Kits returned from Chifeng City (Tested at March 9, 2020)

Specimen

Test results of LYHER® Kit when testing IgG Test results of LYHER® Kit when testing IgM  
 2002003 2002004 2002005 2002003 2002004 2002005

N01 Negative Negative Negative Negative Negative Negative

N02 Negative Negative Negative Negative Negative Negative

N03 Negative Negative Negative Negative Negative Negative

N04 Negative Negative Negative Negative Negative Negative

N05 Negative Negative Negative Negative Negative Negative

P01 Positive Positive Positive Positive Positive Positive

P02 Positive Positive Positive Positive Positive Positive

P03 Positive Positive Positive Positive Positive Positive

P04 Positive Positive Positive Positive Positive Positive

P05 Positive Positive Positive Positive Positive Positive

M-L1 NA NA NA Positive Positive Positive

M-L2 NA NA NA Positive Positive Positive  
M-L3 NA NA NA Positive Positive Positive  
M-L4 NA NA NA Positive Positive Positive  
M-L5 NA NA NA Positive Positive Positive  
G-L1 Positive Positive Positive NA NA NA  
G-L2 Positive Positive Positive NA NA NA  
5

G-L3 Positive Positive Positive NA NA NA  
G-L4 Positive Positive Positive NA NA NA  
G-L5 Positive Positive Positive NA NA NA

The above test results indicate that the transportation has no adverse effect on the product quality of LYHER® Kit.

Table5.2-3: Test results of the expired LYHER® Kit returned from Haikou (will be tested after the expiration date)

Table5.2-4: Test results of the expired LYHER® Kit returned from Chifeng (will be tested after the expiration date)